

Allen.

C-Flex® Head Positioning System



Just Squeeze and Position... Easy Precision Adjustments

The C-Flex head positioner is used for all types of spine procedures, and can be used on an operating table, Jackson™ Table or Allen Flex Frame.

The C-Flex system is the only head positioner on the market that provides the surgeon with full control of the positioning. The surgeon can easily make precision adjustments quickly, just by squeezing the handles and moving the head. In addition, the C-Flex system works with a variety of accessories including both standard and radiolucent Mayfield and DORO skull clamps as well as horseshoes.

Ideal for use with:

- Cervical spine
- Lumbar spine
- Odontoid Screw
- XLIF and lateral procedures
- Fusions
- Laminectomies
- Artificial Disc

At a glance

- Fast precise adjustments of the head
- Cervical flexion / extension
- Head positioning for lumbar surgeries
- Adjust vertically for kyphotic patients
- Works with the Mayfield® Skull Clamp, Mayfield Horseshoe, and DORO® Headrest System
- Includes fold-away mirror
- Easy airway access
- Fits the Allen Flex Frame
- Fits the Jackson Table
- Fits operating tables with Allen C-Flex® AP-Vantage System



A **Hill-Rom** Company

About the Product

The modular system can be used in all types of spine surgery. The quick-connect coupler makes connecting to a Mayfield® skull clamp a snap! Ability to “squeeze and position” supports readjustments intraoperatively. An extensive range of motion makes the C-Flex system ideal for use with extremely kyphotic patients, cervical fusions, and lateral procedures.

Make all of your adjustments from above the table, no need to crouch down or climb under the table as with other positioners on the market.* The fold away mirror enables easy patient monitoring of eyes and airway without having to look under the table.

Technical Specifications

Patient Weight Capacity	227 kg (500 lbs)
Device Weight.....	9 kg (20 lbs)
Lateral Adjustment Range	15 cm (6")
Height Adjustment Range	71 cm (28")
Extension Adjustment Range	23 cm (9")
• +/- 15° lateral tilt	

Allen C-Flex Head Positioning System Contents:

- C-Flex Base Unit
- Prone Mask Head Module
- Flat Plate Head Module
- Starburst Adaptor
- Comfort Mask Disposables
- C-Flex Clean Cape Disposables



Allen C-Flex Base Unit



Prone Mask Head Module



Flat Plate Head Module



Starburst Adaptor



C-Flex Clean Cape Disposables



Comfort Mask Disposables

A-70700 Allen C-Flex Head Positioning System

*Based on comparison to OSI Cervical Management Base Unit, and Mayfield positioners.

Allen® is a registered trademark of Allen Medical Systems. Jackson™ is a trademark of Mizuho Orthopaedic Systems, Inc. Mayfield® is a registered trademark of SM USA, Inc. DORO® is a registered trademark of Pro Med Instruments, Inc. XLIF® is a registered trademark of NuVasive, Inc.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2011 Allen Medical, A Hill-Rom Company.
ALL RIGHTS RESERVED.

D-770591-A2



Sales +44 (0) 870 112 0769
Customer Service +1 978 266 4200 ext. 4286
www.allenmedical.com